

English

1. Introduction

Thank you for having decided for the Illuminated Mini Massager.
Allow yourself a relaxing massage from head to toe.

2. Main features of this device

- Provides a relaxing feeling of well being
- Increases circulation
- Perfect for the Neck, Shoulders, Back, Arms and Legs
- One push of the power button and the three illuminating massage balls start vigorously vibrating
- Ergonomic design
- Can be operated with one hand
- Light weight and handy
- Battery operated (3 x 1,5V AAA)

3. Use

Putting in or changing batteries:

1. Unscrew battery cover by turning it anti-clockwise and remove it.
2. Insert batteries (3x 1,5V Type AAA).
Take care of correct polarity according to picture.
3. Put on battery cover and close it by turning it clockwise.



ATTENTION! – Turn battery cover beyond the “cracking point“. Only then the appliance can be switched on!

4. Switch the appliance on or off by pressing the ON/OFF button.



4. Care and cleaning

Use non-irritant cleaning agents only for cleaning the massager. Wipe off with a damp cloth.

5. Packaging contents

Illuminated Mini Massager
3 batteries 1,5V Type AAA
Instruction manual

6. Technical data:

Model/Type	BA-1056
Voltage:	4,5V
Batteries:	3x 1,5V Type AAA

